







INDUSTRY-SPONSORED SESSIONS

| | | WEDNESDAY, MAY 29 | | | | |
|-------|---------------|-------------------|-------|--|--|--|
| | Congress Hall | South Hall 3AB | | Chamber Hall | | |
| 12:00 | | | | | | |
| | | | 12:15 | Johnson & Johnson Symposium: From evidence to action. SSI reduction strategies | | |
| | | | | Welcome Andrew De Beaux | | |
| | | | | Evidence supporting SSI reduction bundles Marja Boermeester | | |
| | Lunch break | | | Implementation and surveillance of SSI reduction bundles Giles Bond-Smith | | |
| | | | | Industry/clinician partnership to optimise patient outcomes Walter Danker | | |
| | | | 13:15 | Discussions All | | |
| 13:30 | | | | | | |

| | | THURSDAY, MAY 30 | | |
|-------|--|------------------|-------|---|
| | Congress Hall | South Hall 3AB | | Chamber Hall |
| 00 | B. Braun Lecture | | | |
| | Three-Year Follow-up Results of Short Stitch versus Long Stitch Technique in Elective Midline Closure in a Randomized-controlled (ESTOIH) Trial | | | |
| ~~ | René Fortelny | | | |
| 30 | P. Durana I. a strang | Coffee break | | |
| 00 | B. Braun Lecture | | | |
| | Use of prophylactic mesh in contaminated fields | | | |
| | Javier Lopez Monclus | | | |
| | • | | | |
| :30 | | | | |
| 11:45 | Medtronic Symposium: Unlocking potential – optimising the inguinal hernia pathway | | 11:45 | Cousin Surgery Symposium: L is better – from repair to regeneration |
| | Introduction Gill Tierney | | | Moderators: Yohann Renard, Guillaume Pa |
| | Optimising training of surgical trainees Katherine Smith | | | Scientific introduction about PLLA Mehdi Daoudi |
| | Optimising patient throughout Husam Ebied | | | Clinical results of 4DVentral from the Fren Registry (Club Hernie) |
| | Optimising patient care using LPG Yücel Cengiz | Lunch break | | Guillaume Passot (on behalf of David Moszkowicz) |
| | Q&A All | | | Clinical results of 4DMesh from the French Registry (Club Hernie) Yohann Renard |
| | | | | New Italian clinical study on 4DMesh Alberto Sartori, Andrea Sanna |
| | | | | |
| :45 | | | 12:45 | Umbilical hernia repair with 4DUmbilical b PUMP technique Phillip Kirchhoff |